

STERIMAX

2770 Portland Drive, Oakville, Ontario, Canada L6H 6R4

Tel.: 905-890-0661 800-881-3550 • Fax: 905-890-0508 877-546-7667 • Web: www.sterimaxinc.com

July 26, 2022

Compatible Connection Devices for SteriMax Calcium Chloride Injection BP, 100 mg/mL, (DIN 02495945)

Our Contract Manufacturing Organization (CMO) has provided the following list of connection devices that are compatible with the SteriMax Calcium Chloride Injection BP drug product pre-filled syringe (DIN 02495945).

- Baxter Healthcare: Clearlink needle-free luer activated devices
- B. Braun Medical - Ultrasite/Safeflow Connectors
- CODAN Ltd - SWAN-LOCK® Connectors
- ICU Medical - “The Clave connectors that require the use of the 011-CS25 adapter when used with glass syringes include:
 - Clave™
 - MicroClave™
 - MicroClave Clear™
 - NanoClave™
 - Clave™
 - Neutron™
- The ICU Medical connectors that do not incorporate Clave technology and an internal cannula, and are compatible with mating luers or syringes that have an internal diameter range outside of 0.062” to 0.110” (no adapter required) include:
 - Nuitiv™ Clear
 - Tego™
 - Neusite™

DocuSigned by:

Mike Riley

Jul 26, 2022

7336B683D662457...

Mike Riley,
Director, Quality Risk Management & Investigations